


Part Number: **C10PU-MCA2A09***not to scale*

Full Description: Catheter for mouse carotid artery, soft 1Fr-2Fr, 8cm, 9mm insertion length (a0.8), sq tip. LS25, SP25/12.

Unit: each **Pouch Desc:** Mouse CAC PU 1-2Fr 9mm **Pouch:** 3x7.5

Approx Dead Vol: 9 μ L **Connects to:** 25ga **EO Sterilized:** ☒

Intended Species: mouse **Vessel:** carotid artery

Usage Notes:

TUBING

Tip: square **Length 1** **Length 2**

L1: 2.5 cm **L2:** 5.5 cm

Material (IDxOD): PU 1Fr soft (.18x.35 $\pm$.05mm) **Material (IDxOD):** PU 2Fr (.43x.69 $\pm$.03mm)

FEATURES

Distal (tip) end

Feature	Description	Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 0.8 cm	☒
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 2.5 cm	☒
Feature C:		c: cm	☐

Proximal end

Feature	Description	Location	Fixed
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1: LS25 **Part 2:** SP25/12