


Part Number: **C10PU-MCA2A10***not to scale*

Full Description:	Catheter for mouse carotid artery, soft 1Fr-2Fr, 8cm, 10mm insertion length (a0.9), sq tip. LS25, SP25/12.		
Unit:	each	Pouch Desc:	Mouse CAC PU 1-2Fr 10mm
		Pouch:	3x7.5
Approx Dead Vol:	9 μ L	Connects to:	25ga
		EO Sterilized:	☒
Intended Species:	mouse	Vessel:	carotid artery
Usage Notes:			

TUBING

Tip:	square	Length 1	Length 2
		L1:	L2:
		2.5 cm	5.5 cm
		Material (IDxOD):	Material (IDxOD):
		PU 1Fr soft (.18x.35 $\pm$.05mm)	PU 2Fr (.43x.69 $\pm$.03mm)

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 0.9 cm	☒
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 2.5 cm	☒
Feature C:		c: cm	☐
Proximal end			
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1:	LS25	Part 2:	SP25/12
---------	------	---------	---------